

MOBILE COMMODE - SPLIT SEAT - STAINLESS STEEL

freedomhca
HEALTHCARE

www.freedomhc.co

Quality System
Quality
Endorsed
Company
ISO 9001 : 2008

- The new seating provides full access to all the body parts of the patient.
- Absolutely no need for the user to stand off the chair to have their body parts washed.
- The frame is made from 304 marine grade stainless steel.
- Life time warranty against corrosion.
- Made tough. Expect a service life of around 10 years minimum.
- Suitable for rental or domestic use.
- For hospitals, nursing homes and other institutional use.
- Low maintenance.

PATENTED

HBA355

New ARMREST



- multiple height adjustable
- lowers to the seat level for easy patient transfer
- padded armrest

- Backrest •
- contoured backrest
 - rubber padded for comfort
 - rubber surface prevents coldness upon skin contact

LOOK!



SEAT Patented - New



- the new seating allows the patient to hygienically wash their body parts without getting off the chair.
- The front and the rear of the seat has an open access cut out.
- it allows the carer to clean the body parts of the patient from the front and the rear.

Swivel armrest

- Foot rest •
- slides forward and backward
 - rubber padded for comfort
 - rubber surface prevents coldness upon skin contact

- Frame •
- stainless steel marine grade 304

Optional BACKREST



- TSA350 - L Reclining backrest.
- approx. 20 degree
 - suitable for washing hair

New HAND GRIP



The handle grips rotate to secure the patient and prevents the patient from falling forwards.

- Castors
- all lockable
- low maintenance

Patented Feature

The patient does not need to stand up, off the chair



Rear Front
Full access to body parts Full access to body parts

W

1 YEAR WARRANTY
* conditions apply

HBA355

overall height	overall length	overall width	seat depth	seating width	seat height	armrest height	footrest length	castor size	unit weight	weight capacity
104cm	67cm	57cm	42cm	48cm	58cm	min0cm max31cm	min0cm max26cm	12.5m	17.9kg	200kg

Quality Assurance

Tested by
a NATA
accredited
laboratory

AS/NZS
3973 : 2009